


ATLANTIS™

ATLANTIS™ abutments












– for all major implant systems

Material
reference guide



DENTSPLY
IMPLANTS

Biocompatible materials, ATLANTIS™ abutments

	Titanium	Titanium gold-shaded	Zirconia				
	Milled from solid titanium alloy blanks, titanium Alloy 6Al-4V (Grade 5)		Milled from solid zirconia blanks. Available in five shades, including white and translucent white. Yttria-stabilized tetragonal zirconiapolycrystals (Y-TZP)				
		Thin coating of titanium nitride (TiN) which gives a warm gold-shaded color.	Shade 00T has a higher translucency allowing for more light transmission.				
			shade 00T	shade 00	shade 10	shade 20	shade 30
Guideline for usage			Used with very light and bleached shades.	Used with very light and bleached shades.	Used with light to mid range shades, related to the Vita shade guide, e.g. B1	Used with mid range shades, related to the Vita shade guide, e.g. A2.	Used with stronger/darker shades, related to the Vita shade guide, e.g. B3.
ATLANTIS™ Abutment							
for cement-retained restorations			N/A				
ATLANTIS™ Crown Abutment							
for single-tooth, screw-retained restorations	N/A	N/A					

NOTE: Images shown for communication purposes only and should not be used as a representation of the exact material color.